

# PREMIUM 400

## 4 CHANNELS ELECTROTHERAPY | 258 PROGRAMS

OVER 250 PROGRAMS IN YOUR HANDS, PREMIUM 400 IS THE IDEAL ELECTROSTIMULATOR FOR THE SPORT AND THE CARE OF THE BEAUTY, FOR MEN AND WOMEN.

58 FITNESS | 53 SPORT | 44 SPECIALS SPORT | 36 BEAUTY | 23 MICROCURRENTS | 18 SERIAL SEQUENTIAL STIMULATION  
12 TENS AND PAIN | 7 ACTION NOW | 3 INCONTINENCE | 1 IONOPHORESIS

### TECHNICAL FEATURES

#### Display

Visible area size 2,6"

#### Channels

4 independents (8 electrodes)

#### Frequency

0,3-150Hz \*

#### Pulse amplitude

40-450µs \*

#### Power

0-120mA per channel

#### Power supply

Rechargeable batteries

#### Size

mm 100x160x35


#### Weight

gr 454

\* According to the type of current

### EQUIPMENT

- 1 Bag
- 1 Stimulator PREMIUM 400
- 4 Cables for electrodes connection
- 2 Cables for microcurrents and ionophoresis
- 4 Self-adhesive square electrodes
- 4 Self-adhesive rectangular electrodes
- 1 Charger
- 1 Operating manual



CE  
0476



TENS



MENS



IONO



3S



ACTION NOW



RUN TIME



LAST 10



FAVORITES



MULTI USER



PROGRAMMABLE



STIM LOCK



FUNCTION 2+2



MY TRAINER



SYNCHRO STIM



WORK TIME



RECHARGEABLE BATTERIES

# PREMIUM 400

## 4 CHANNELS ELECTROTHERAPY | 258 PROGRAMS

### PROGRAMS LIST

#### SPORT

Capillarization

Warm-up

Pre-competition warm-up

Active recovery

Maximum strength

Endurance strength

Explosive strength

Aerobic endurance

Reactivity

Post-competition recovery

Decontracting

Hypertrophy

#### SPECIAL SPORT

Running

Cross-country skiing

Martial arts

Tennis

Soccer

Cycling

#### FITNESS

Firming

Bio-Pulse firming

Sculpting

Bio-Pulse sculpting

Toning

Mass building

Body sculpting

Definition

Jogging

Anaerobic fitness

Aerobic fitness

Cramp prevention

#### BEAUTY

Drainage

Bio-Pulse drainage

Lipolysis

Post-pregnancy lipolysis

Toning massage

Connective massage

Swollen arms

Face capillaries

Skin tone improvement

Post-pregnancy drainage

Post-pregnancy firming

Breast firming

Breast sculpting

Lifting effect

#### REHAB

Quadriceps atrophy

Recovery after ACL surgery

Shoulder subluxation prevention

#### TENS and PAIN

Conventional antalgic tens

Modulated antalgic tens

Endorphinic tens

Menstrual pain

Chronic pain

Cervical pain

Muscle pain

Knee pain

Scapulohumeral periarthritis

Chronic lumbago

Bursitis-tendinitis

Osteoarthritis

#### MICROCURRENTS

Epicondylitis

Periarthritis

Muscle restoration

Contusion

Hedema

Skin ulcera

Sciatica

Lumbago

Brachial neuralgia

Acute pain

Articular pain

Stiff neck

Whiplash

Cervical spondylosis

Shoulder sprain

Knee sprain

Carpal tunnel

Osteoarthritis

Ankle Sprain

Achille tendon inflammation

Patella tendon inflammation

Rotator cuff inflammation

Tendon Inflammation

#### INCONTINENCE

Mixed incontinence

Stress incontinence

URGE incontinence

#### IONOPHORESIS

Ionophoresis

#### SERIAL SEQUENTIAL STIMULATION

The "3S" program list includes 18 parameter combinations.

The "3S" programs are characterized by a delay activation of the channels 3 and 4 compared with the channels 1 and 2. The Serial Sequential Stimulation permits to stimulate the musculature in kinetic chain thanks to the differentiated activation times of the muscular groups involved.

#### ACTION NOW

The Action Now program list includes 7 parameter combinations.

This mode is suggested in sport field for athletic preparation where you wish to add the muscular contraction induced by a stimulator to a work made with overlads in dynamic and isometric form.