# **MAGNUM 3000 PRO**

# PROFESSIONAL MAGNETOTHERAPY | 70 PROGRAMS

POWERFUL PROFESSIONAL MAGNETOTHERAPY EQUIPPED WITH MANY SPECIAL PROGRAMS AND FUNCTIONS THAT MAKE IT EFFECTIVE THE TREATMENT OF SEVERAL OSTEOARTICULAR PATHOLOGIES AND THE RICOVERY AFTER ORTHOPEDIC SURGERY. POSSIBILITY OF LOCAL AND TOTAL-BODY TREATMENTS.

## **TECHNICAL SPECIFICATIONS**

## Display

Size of the visible area 2.6"

#### Channels

2

#### Intensity can be emitted with flexible diffuser

- 200 Gauss on average for channel
- 400 Gauss on average on the whole (for 2 channels)

#### Frequency

5-200 Hz

#### **Duty Cycle**

5-50%

#### Free programming

Possibility of adding up to 10 custom programs

#### 1+1 function

Possibility of running 2 different programs at the same time, using 1 channel each

#### Power supply

- Mains
- Rechargeable batteries (not included)

#### Sizde | Weight

170x220x60 mm | 820 gr








## STANDARD EQUIPMENT

#### Code **G5335**

- · 1 Bag
- · 1 MAGNUM 3000 PRO magnetotherapy device
- · 2 "Flexible" solenoids, 300x100 mm
- · 2 Elastic bands with Velcro 800x80 mm
- · 1 Magnum Test
- · 1 Power supply unit
- · 1 Operating manual













## **PURCHASE OPTION**

Code **G6034** MAGNUM 3000 PRO with 2 solenoids "Pocket Pro" and 1 elastic band with Velcro



EASY PROGRAM



1+1 MODE



LAST 10



PROGRAMMABLE



MULTI USER



**RENTAL OPTION** 



**RUN TIME / EDIT PARAMETERS** 



MAINS SUPPLY



# **MAGNUM 3000 PRO**

# PROFESSIONAL MAGNETOTHERAPY | 70 PROGRAMS

# **PROGRAM LIST**

ARTHROSIS/ARTHRITIS	NIGHT PROGRAMS
Knee arthrosis	Night program 15 Hz
Cervical arthrosis	Night program 50 Hz
Arthritis	Night program 75 Hz
Arthrosis	DEFAULT PROGRAMS WITH DIFFERENT FREQUENCY
Hand arthrosis	Multi-frequency program
Hip arthrosis	Program 5 Hz
Shoulder arthrosis	Program 10 Hz
Column arthrosis	Program 15 Hz
PAIN	Program 20 Hz
Algodystrophy - Sudek disease	Program 25 Hz
Back pain	Program 30 Hz
Joint pain	Program 35 Hz
INFLAMMATIONS	Program 40 Hz
Anti-inflammatory	Program 45 Hz
Tendinitis	Program 50 Hz
Epicondylitis	Program 55 Hz
Epitrochleytis	Program 60 Hz
Carpal tunnel	Program 65 Hz
Shoulder tendinitis	Program 70 Hz
ORTHOPEDICS	Program 75 Hz
Post ACL surgery	Program 80 Hz
Cartilage lesions	Program 85 Hz
Knee posthesis	Program 90 Hz
Hip prosthesis	Program 95 Hz
Spinal fusion	Program 100 Hz
Osteotomy	Program 120 Hz
BONE PATHOLOGIES	Program 140 Hz
Fractures	Program 160 Hz
Scaphoide fracture	Program 180 Hz
Wrist fracture	Program 200 Hz
Shoulder fracture	<u></u>
Collarbone fracture	
Forearm fracture	<u></u>
Pelvic fracture	
Tibial fracture	
Femoral fracture	<u></u>
Delayed fracture consolidation	<u></u>
Fractture with plaster	



TRAUMAS / EDEMAS / WOUNDS

Osteonecrosis
Osteoporosis
Pseudoarthrosis

Ankle sprain

Edemas- contusions

Ulcers

