

DIACARE 5000 | DIACARE 5000 RE

TECARTHERAPY | 14 PROGRAMS

DIACARE 5000 INCLUDES ALL THE GLOBUS TECARTHERAPY QUALITY IN A DEVICE WITH REDUCED SIZE AND WEIGHT, TO GUARANTEE THE MAXIMUM PORTABILITY, MAINTAINING THE SAME EFFECTIVENESS AND PROFESSIONALISM.

TECHNICAL FEATURES

Display

Size 4.3" color touch screen display

Nominal Power (NP)

180W ±10%

Emission frequency

300kHz, 400kHz, 448kHz, 470kHz, 700kHz, 1000kHz, ± 10%

Treatment modality

Capacitive and resistive

Capacitive electrodes

3 size Ø (mm 30, 50, 70)

Resistive electrodes

3 size Ø (mm 30, 50, 70)

Mains supply

100-240V AC, 50-60Hz, 1,8A

Size

170x220x60 mm

Weight

1,04 Kg

EQUIPMENT

· 1 Carrying case

· 1 DIACARE 5000 tecartherapy unit

· 1 Counter electrode (plate)


· 1 Monopolar resistive handpiece with 3 heads
Ø (mm 30, 50, 70)

· 1 Monopolar capacitive handpiece with 3 heads
Ø (mm 30, 50, 70)

· 1 Cream can ml 250

· 1 Power supply unit

· 1 Operating manual



PURCHASE OPTIONS

· Cod. **G3488** DIACARE 5000

· Cod. **G3996** DIACARE 5000 RE *

* DIATHERMIE RECHARGEABLE:

- 5 HOURS OF TREATMENT INCLUDED FREE OF CHARGE
- AFFORDABLE RE-CHARGING COST
- GUARANTEED ASSISTANCE FOR 5 YEARS AND BACKUP EQUIPMENT
- RECHARGE ON WEBSITE



DIACARE 5000 | DIACARE 5000 RE

TECARTHERAPY | 14 PROGRAMS

PROGRAMS LIST

MEDICAL PROGRAMS

Muscle trauma pain

Tendon trauma pain

Ankle traumas

Lumbagos

Shoulder pain

Cervical pain syndromes

Pulled muscles

Joint diseases

Arthrosis of the knee joint

Whiplash

Sciatica

Slipped disc

Piriformis syndrome

Postoperative treatment of femur fractures



PRESET PROGRAMS



PULSED TREATMENTS



ATHERMAL TREATMENTS



MULTI-FREQUENCY



AUTOMATIC FREQUENCY SCANNING



MYOFASCIALS



HANDS-FREE MODE



TRANSDERMAL VEHICULATION



RESISTIVE



CAPACITIVE



BIPOLAR



COLOR DISPLAY



TOUCH SCREEN DISPLAY



MAINS SUPPLY



GLOBUS RECHARGE

